


Investigation of Hits against the GLP-1 Receptor Identified by Affinity Selection Mass Spectrometry: Hit Expansion, Structure-Activity Relationships, and Selectivity Profiling

Zane Thistleford¹, Christina Thompson¹, William LaMarr¹, Can Özbal¹, Jan Kubicek², Barbara Maertens²

¹Momentum Biotechnologies, 3 Federal St., Suite 300, Billerica, MA 01821 ²Cube Biotech, Creative Campus Allee 12, 40789 Monheim, Germany



momentum
BIOTECHNOLOGIES

Abstract


A prior high-throughput affinity selection mass spectrometry (ASMS) screen of 50,000 compounds against the GLP-1 receptor expressed in copolymer-based nanodiscs using the NativeMP™ platform by Cube Biotech identified 34 confirmed hit compounds. Here we report the hit expansion and characterization of three distinct chemical series emerging from that screen. Structurally similar analogs were acquired from commercial sources. Three chemical scaffolds were characterized: a pyrazolinone/indazolone lactam series (20 analogs), an aminopyrimidine series (16 analogs), and a benzimidazole-amide series (11 analogs). Key structure-activity findings included: a 5–100x potency advantage for direct C4-aryl over aryl-ketone substituents in the lactam series; a 26-fold regiochemistry penalty for meta- vs. para-substitution on the pendant phenyl; dominant halogenation effects in the aminopyrimidine series ($\text{Br} \approx \text{CF}_3 > \text{Cl} \gg \text{H}$); and a 3x activity gain from a constrained piperidine linker in the benzimidazole-amide series. Selective GLP1R binders were subsequently profiled against two related off-targets GIPR and GCGR all in nanodisc format prepared by Cube Biotech. GCGR was the dominant cross-reactivity challenge (39% of compounds preferred GCGR). Halogenated compounds showed 3x better selectivity on average, and competition assays confirmed orthosteric binding within each series. This work illustrates the utility of nanodisc-based ASMS for rapid, label-free hit expansion and selectivity differentiation of membrane protein binders.

Background & Objectives


A previously described and presented ASMS screen of the GLP-1 Receptor identified 34 confirmed hits from a 50K compound diversity library. This study builds on that work by:

- Hit expansion of select chemotypes via chemical clustering (Butina, Tanimoto ≥ 0.20)
- SAR analysis using ASMS peak area as a proxy for binding affinity
- Selectivity profiling vs GIPR, GCGR and blank nanodisc


GLP1R: Type 2 Diabetes and Obesity Target



ALIS-ASMS Workflow



Affinity Determination — Dose-Response Curves



Example concentration response curves for confirmed hit from the ASMS HTS campaign comprising 50,000 diversity compounds from Momentum P57 chemical library


Compound Clustering Summary

	N	Hit	Scaffold	Area Range
A	20	1271517	Pyrazolinone/Indazolone Lactam	0.4K–103K
B	16	1038635	Aminopyrimidine (chloro-aryl)	10K–2.1M
C	11	1293841	Benzimidazole-amide	0.4K–104K

Three HTS Hits Selected as Starting Points for Hit Expansion

★ 1271517 20 analogs A — Pyrazolinone series Pyrazolinone/Indazolone lactam Area: 42.5K Score: 98.4	★ 1038635 16 analogs B — Aminopyrimidine series Aminopyrimidine (chloro-aryl) Area: 190.5K Score: 99.3	★ 1293841 11 analogs C — Benzimidazole series Benzimidazole-amide Area: 33.8K Score: 98.7
---	--	---

Cluster A — Pyrazolinone / Indazolone Lactam Series (20 Cpd)



★ 1271517 ASMS Area: 42.5K 1536202	1151192 ASMS Area: 103.1K 1273111	1266629 ASMS Area: 79.4K 1662436	1273094 ASMS Area: 78.7K 1110901
--	---	--	--

 ASMS Area: 69.9K 1662434	 ASMS Area: 55.6K 1537119	 ASMS Area: 40.8K 1537177	 ASMS Area: 40.4K 1537131
--	--	--	--

 ASMS Area: 35.3K 1541838	 ASMS Area: 30.9K 1000977	 ASMS Area: 25.2K 1536994	 ASMS Area: 24.5K 1537105
--	--	--	--

 ASMS Area: 7.1K 1537169	 ASMS Area: 3.6K 1271520	 ASMS Area: 2.8K 1537138	 ASMS Area: 2.2K 1266624
---------------------------------------	---------------------------------------	---------------------------------------	---------------------------------------

 ASMS Area: 1.9K	 ASMS Area: 1.6K	 ASMS Area: 415	 Non-binder
---------------------	---------------------	--------------------	----------------

Clusters B & C — Hit Expansion Overview

Cluster B — Aminopyrimidine Series (16 compounds, incl. 1 non-binder)
16 analogs screened around original hit 1038635. Halogenation is the dominant SAR driver — para-Br and CF₃ deliver 10–50x improvement. Best compound: 1104065 (diphenylamino, 2.1M area, best in class). Regiochemistry critical: ortho-OMe/para-Cl (514K) vs swapped isomer (42K) = 12x penalty.

Cluster C — Benzimidazole-Amide Series (11 compounds, incl. 2 non-binders)
11 analogs screened around original hit 1293841. Best compound: 1272018 (piperidine-constrained linker, 104K, 3x gain over hit). Benzimidazole core irreplaceable (pyrrole swap = 80x loss). EWG at C5/C6 required; OMe → 5x loss. 2 non-binders confirm scaffold essentials.

Structure–Activity Relationships: Cluster A (Pyrazolinone/Indazolone Lactam)

Two Distinct Sub-Series Within Cluster A
Sub-series 1: Direct aryl at C4 — top binders (3.6K–103K). **Sub-series 2:** Aryl ketone at C3 — weaker (0.4–31K). Direct-aryl is 5–100x more potent.

★ Top binder 1273094 Area: 78.7K Score: 99.6	Top binder 1536202 Area: 69.9K Score: 97.6	★ Original Hit 1271517 Area: 42.5K Score: 98.4
Ketone sub-series 1537119 Area: 30.9K Score: 98.4	Meta-Cl (weak) 1271520 Area: 1.6K Score: 72.0	Weakest binder 1537138 Area: 415 Score: 10.3
C4 direct aryl essential — ketone penalised 5–100x Direct C–C bond to pendant aryl at C4 is 5–100x more potent than a ketone linker. Within the ketone series, ortho-substitution rescues activity: ortho-OMe (31K) vs unsubstituted (2.8K) via preferred rotamer.	Para-substitution on pendant phenyl strongly preferred Para-OAllyl (79K) and para-Cl are optimal. Meta-Cl (1271520) drops 42K→1.6K — a 26x loss from a single regiochemical change. Pocket is stringently defined.	Para-Cl on indazole ring is consistently additive (~2x) Adding para-Cl to the indazole provides ~2x gain regardless of pendant aryl substitution. Likely a discrete halogen-bond or lipophilic contact in a sub-pocket.
Ortho-substitution rescues ketone sub-series Ortho-OMe ketone (1537119, 31K) is 11x better than unsubstituted phenyl ketone (1536994, 2.8K). Ortho group enforces a dihedral mimicking the direct-aryl geometry.	High TPSA penalised — hydrophobic binding pocket Top binders: avg TPSA 57 Å² vs 76 Å² bottom quotient. H-bond donors also lower (1.5 vs 2.3). Consistent with a hydrophobic, solvent-excluded binding site.	

Structure–Activity Relationships: Cluster B (Aminopyrimidine Series)

Best in class 1104065 Area: 2.11M 99.8	CF₃ (top 3) 1006677 Area: 69.19K 98.6	Para-Br 1006653 Area: 537.0K 99.5
Correct regio. 1024731 Area: 514.8K 99.2	Swapped regio. ⚠ 1006673 Area: 42.0K 96.7	No halogen 1235165 Area: 32.6K 98.6
Halogenation is the dominant driver Br = CF ₃ > Cl > H. 10–50x gains over unhalogenated. Best compound 1104065 (diphenylamino, 2.1M) is the highest area in the entire dataset.	Regiochemistry: 12x penalty for a simple swap Ortho-OMe/para-Cl (1024731, 514K) vs para-OMe/ortho-Cl (1006673, 42K). Same atoms, different positions — binding pocket is spatially stringent.	Bicyclic tail sub-series: distinct chemotype 1095779, 1147707, 1080458 lack the aniline-pyrimidine pharmacophore yet bind well — may access the same site via an alternative mode.


Structure–Activity Relationships: Cluster C (Benzimidazole-Amide Series)

Best in class 1272018 103.8K 99.8	★ Original Hit 1293841 33.8K 98.7	F analog 1335417 31.7K 98.3
OMe (weak) ⚠ 1272495 7.1K 75.5	Pyrrole (dead) 1253450 432 14.3	Lead — constrained linker 1272018b 103.8K 99.8
Constrained linker delivers 3x gain 1272018 (103K) exceeds hit 1293841 (34K) 3x. Piperidine ring locks N-methyl indazole in optimal geometry — highest-priority lead	EWG at C5/C6 required; EDG kills activity Cl and F are equipotent at C5. OMe at C6 (1272495, 7K) is 5x weaker — electron donation reduces N-mediated binding contacts.	Benzimidazole core is irreplaceable Pyrrole swap (1253450) = 432 area — ~80x collapse. Fused bicyclic aromatic system essential for π-stacking or H-bond acceptor contacts

GLP1R Selectivity Screening: GLP1R, GIPR, GCGR

GLP1R, GIPR, and GCGR are three closely related Class B1 GPCRs sharing a common architecture: a large N-terminal extracellular domain (ECD) for peptide hormone recognition, a seven-transmembrane (7-TM) helical bundle, and an intracellular tail mediating Gs protein coupling. The receptors share substantial transmembrane sequence identity: approximately 44% between GLP1R and GIPR, and ~48% between GLP1R and GCGR reflecting their common evolutionary origin. The high structural similarity of their transmembrane binding pockets, especially between GLP1R and GCGR is the primary reason selective small molecule engagement of any one receptor over the others remains a central challenge in drug discovery for this target family.

Target Receptor Structures — 7-TM Topology (GLP1R, GIPR, GCGR)



45 compounds screened in parallel as singletons against GLP1R, GIPR, GCGR and blank nanodisc. Selectivity = log₂(GLP1R/off-target). Values >0 = GLP1R preferred. **22/44 compounds GLP1R-selective. GCGR is the dominant challenge (39% of compounds prefer GCGR).**

Top GLP1R-Selective Compounds

1311531 composite: 12.30 vs GIPR: 0.8, vs GCGR: 18.1, vs Blank: 18.1 LogP 4.28, TPSA 46.1, MW 321.4, cFSP3 0.44, Halo 1 GLP1R 275K, GIPR 163K, GCGR 0, Blank 0	1299341 composite: 5.72 vs GIPR: 0.5, vs GCGR: -2.3, vs Blank: 18.9 LogP 3.45, TPSA 65, MW 389.9, cFSP3 0.35, Halo 1 GLP1R 499K, GIPR 354K, GCGR 2.38M, Blank 0	1271520 composite: 3.23 vs GIPR: 0.6, vs GCGR: 2.2, vs Blank: 6.9 LogP 4.2, TPSA 57.8, MW 323.8, cFSP3 0.11, Halo 1 GLP1R 578K, GIPR 385K, GCGR 125K, Blank 5K	1660920 composite: 3.19 vs GIPR: 1.0, vs GCGR: 2.8, vs Blank: 5.8 LogP 3.86, TPSA 57.8, MW 303.4, cFSP3 0.15, Halo 0 GLP1R 499K, GIPR 256K, GCGR 73K, Blank 9K
1293685 composite: 3.08 vs GIPR: 0.1, vs GCGR: 1.5, vs Blank: 18.1 LogP 3.9, TPSA 47.6, MW 369.4, cFSP3 0.32, Halo 1 GLP1R 268K, GIPR 197K, GCGR 33K, Blank 5K	1273451 composite: 3.00 vs GIPR: 0.3, vs GCGR: -3.6, vs Blank: 0.3 LogP 3.13, TPSA 49.8, MW 369.4, cFSP3 0.47, Halo 0 GLP1R 241K, GIPR 110K, GCGR 2K, Blank 106K	1663759 composite: 2.86 vs GIPR: 0.5, vs GCGR: 1.5, vs Blank: 6.9 LogP 5.09, TPSA 66.3, MW 433.5, cFSP3 0.35, Halo 1 GLP1R 219M, GIPR 1.54M, GCGR 794K, Blank 22K	1275764 composite: 2.86 vs GIPR: 0.1, vs GCGR: 0.9, vs Blank: 7.2 LogP 3.09, TPSA 41.5, MW 288.4, cFSP3 0.29, Halo 0 GLP1R 274K, GIPR 192K, GCGR 146K, Blank 2K

Least GLP1R-Selective Compounds: High Off-Target Binding

1104701 composite: -3.05 vs GIPR: -4.0, vs GCGR: -5.5, vs Blank: 0.3 LogP 4.04, TPSA 45.3, MW 326.4, Halo 0 GLP1R 2K, GIPR 31K, GCGR 86K, Blank 2K	1299453 composite: -1.57 vs GIPR: 0.3, vs GCGR: -3.6, vs Blank: -1.5 LogP 3.89, TPSA 50.5, MW 391.5, Halo 0 GLP1R 163K, GIPR 133K, GCGR 1.92M, Blank 488K	1298002 composite: -1.22 vs GIPR: -0.0, vs GCGR: -2.4, vs Blank: -1.2 LogP 4.86, TPSA 75.6, MW 373.5, Halo 0 GLP1R 39K, GIPR 40K, GCGR 200K, Blank 89K	1306217 composite: -1.22 vs GIPR: -0.1, vs GCGR: -2.8, vs Blank: -1.4 LogP 4.86, TPSA 64.2, MW 442.5, Halo 1 GLP1R 1.80M, GIPR 1.80M, GCGR 8.83M, Blank 4.02M
--	---	--	---

1. Blank nanodisc is the most informative discriminator

25/44 compounds (57%) bind GLP1R at least 4x more than blank nanodisc. This is the cleanest selectivity signal. Compounds with low blank binding likely bind via a specific, defined mode rather than non-specific lipophilic partitioning into the nanodisc membrane. This axis shows essentially no correlation with any single descriptor, suggesting it is driven by specific binding contacts rather than bulk physicochemical properties.

2. GCGR is the dominant off-target challenge

17 of 44 compounds (39%) show higher GCGR signal than GLP1R (log₂ < 0). GCGR (glucagon receptor) is the closest structural relative of GLP1R, making cross-reactivity unsurprising. LogP shows the strongest (modest) correlation vs GCGR (r=+0.17): higher LogP slightly favours GLP1R selectivity. GCGR-selective compounds have notably lower aromatic ring count (avg 2.1 vs 2.7), suggesting flat, aromatic scaffolds may be less discriminating between these two receptors.

3. Halogen presence is consistently associated with GLP1R selectivity

Halogenated compounds show mean composite selectivity of +1.93 vs +0.63 for non-halogenated (3x difference). The effect is most pronounced in the blank nanodisc axis (avg +5.1 vs +2.5 log₂), suggesting halogens promote specific binding rather than simply increasing lipophilicity. This is consistent with halogen bonding or tight hydrophobic contacts at defined GLP1R binding pockets not present on GCGR/GIPR.

4. Low TPSA / low aromatic ring count favours selectivity

Compounds with TPSA < 50 Å² show mean composite selectivity of +1.81 vs +0.84 for higher TPSA compounds. The negative correlation of aromatic ring count (r = -0.19) is the strongest single descriptor effect: highly aromatic, flat compounds tend to bind non-selectively. This may reflect pan-assay interference or non-specific hydrophobic binding. Compounds with higher sp³ character (cFSP3 ≥ 0.35) average +1.34 vs +0.89 selectivity — 3D character modestly benefits selectivity.

Conclusions

Small molecule engagement of the GLP-1 receptor using affinity selection mass spectrometry (ASMS) proved to be a productive starting point for drug discovery, with a 50,000 compound diversity screen identifying three structurally distinct hit series that could be expanded and characterized in depth. Hit expansion across 52 analogs established clear structure-activity relationships within 3 scaffolds. Selectivity profiling against the closely related GIPR and GCGR confirmed that GLP1R-selective binding is achievable (22 of 44 compounds showed net GLP1R preference) but that GCGR cross-reactivity, driven by ~48% transmembrane sequence identity, represents the primary differentiability challenge. Taken together, this work demonstrates that nanodisc-based ASMS is a powerful and efficient platform for the iterative characterization of membrane protein binders, from primary hit identification through analog profiling and selectivity differentiation, without the need for labeled compounds, surface immobilization, or extensive assay development, and that it is directly applicable to other GPCRs and membrane-embedded targets where selectivity among structurally homologous family members is a key challenge.